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[1] Pantev, C., Okamoto, H., & Teismann, H. (2012). Music-induced cortical plasticity and lateral inhibition in the human auditory cortex as foundations for tonal Tinnitus treatment. *Frontiers in Systems Neuroscience*, 6. DOI: 10.3389/fnsys.2012.00050.

[2] Strauss, D., Corona-Strauss, F., Seidler, H., Haab, L., & Hannemann, R. (2016), Notched environmental sounds: A new hearing aid-supported tinnitus treatment evaluated in 20 patients. *Clinical Otolaryngology*. DOI: 10.1111/coa.12575.

[3] Powers, L., & Munhóes dos Santos, G. (2016). Notch Therapy: a new approach to tinnitus treatment. In submission.

[4] Powers L, Munhóes dos Santos G. (2015). Acoustic Therapy for the Treatment of Tinnitus: A Primer. *Hearing Review*, 22(6). 22.

- [5] Powers, L., & Ramirez, T. (2014). Clinical management of tinnitus. *AudiologyOnline*, Article 12558. Retrieved from: <http://www.audiologyonline.com/articles/clinical-management-of-tinnitus-12558>
- [6] Haab, L., Lehser, C., Corona-Strauss, F., Bernarding, C., Seidler, H., Strauss, D., Hannemann, R. (2016). Six-Month Evaluation of a Hearing Aid Supported Tinnitus Treatment Using Notched Environmental Sounds. In submission.
- [7] Henry, J. Flick, C., Gilbert, A., Ellingson, R., Fausti, S. (2001). Comparison of two computer-automated procedures for Tinnitus pitch matching. *Journal of Rehabilitation Research and Development*, 38(5). 557–566.
- [8] Henry, J. & Meikle, M. (2000). Psychoacoustic measures of tinnitus. *Journal of the American Academy of Audiology*, 11. 138-155.
- [9] Dillon, H. (2012). *Hearing Aids* (2nd ed.). New York: Thieme.
- [10] Mueller, G., Bentler, R., & Ricketts, T. (2014). *Modern Hearing Aids: Pre-Fitting Testing and Selection Considerations*. San Diego: Plural.
- [11] *Progressive Tinnitus Management (PTM) – Level 2*. Veterans' Administration National Center for Rehabilitative Auditory Research (NCRAR).
http://www.ncrar.research.va.gov/Education/Documents/TinnitusDocuments/05_HenryPTMHB_39-50.pdf
- [12] *Tinnitus Functional Index*. Veterans' Administration National Center for Rehabilitative Auditory Research (NCRAR).
<http://www.ncrar.research.va.gov/Education/Documents/TinnitusDocuments/TFI.pdf>
- [13] *Tinnitus Hearing Survey*. Veterans' Administration National Center for Rehabilitative Auditory Research (NCRAR).
<http://www.ncrar.research.va.gov/Education/Documents/TinnitusDocuments/THS.pdf>
- [14] *Self-Efficacy for Managing Reactions to Tinnitus (SMRT)*. Veterans' Administration National Center for Rehabilitative Auditory Research (NCRAR).
<http://www.ncrar.research.va.gov/Education/Documents/TinnitusDocuments/SMRT.pdf>